(English)


PyraBack® was designed specially to:

- 1) Treat back pain , and bad posture
- 2) After several sessions helps you to get back flexibility.
- 3) It can be used while lay down or when sitting.
- 4) It is good for any person that spend several hours sitting.
- 5) It can be used for those that had bad posture or complain of back pain or have limited back movement.
- 6) It will help you to eliminate your continuous back pain
- 7) Provides you with a digit-massage without having to have the inconvenience to move or wait for a specialist appointment.
- 8) It is very easy to carry or use.
- 9) It will alleviate your pain in minutes without having to take medications.

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
Email: Info@pyraback.com

www.PyraBack.com

WHAT IS THE
PYRABACK? /
“¿QUÉ ES EL
pyraback®?”

Treats your Back Pain

**USA / International Patented
& FDA Approved**

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In 1971, McKenzie introduced a new method to treat back pain based in back extension. It took USA's practitioners years to start using this technique and after many years of research they proved its benefits. Actually, the McKenzie treatment is used in almost all cases of back pain. The Pyra-Back® is based not only in the McKenzie Therapy, but the Trigger point theories. This other method of treatment was introduced by Janet Travell, MD. (who was President JFK physician) and David Simmons, M.D. pioneers of "Trigger Point" in print in 1942.

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./.

The use of the Pyra-Back® is very simple, complemented with specific stretching back-exercise. It is placed in a flat surface with the knots

PyraBack®
Superior View /
Vista Superior

up while the person from sitting position lay down allowing the knots enter in contact with his/her Para-spinals muscles, starting on the Lumbar-Costal area. Use it as mush as need it and for a period of no longer of 15 min. Not recommended to be used over naked skin. The above picture shows the most common trigger points in the back.

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MCKENZIE THEORY

The intradiscal pressure decreases as it moves into back extension. Beside the sitting, the lost of back extension and the frequency of flexion produce back pain. Repeated and sustained flexions stresses the nucleous and force it to move posteriorly, and if the disk loose its hydrostatic mechanism could cause a herniation of the nucleous pulposus with nerve root compression. The back extension, reduce the stress on the posterior annular wall and simultaneously cause the fluid nucleous to move anteriorly. ./.

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TRIGGER POINT THEORY

Trigger Point is "a highly irritable localized spot of exquisite tenderness in a nodule in a palpable taut band of muscle tissue." In other words, "nodule" or trigger point hurts when you push on it. The place where the contraction occur is called sarcomeres which act as a tiny pumps contracting and releasing to circulate blood. A knot is a mass of sarcomeres in the state of maximum contraction formed in the muscular fiber. When these lack oxygen and accumulate waste products of metabolism irritate the trigger point, then pain is the result. The treatment is deep pressure-stroking massage. Many Massagist use a finger knot, The knobble massage tool or Thera-Cane®. The PyraBack uses the Knot too and combine it with a passive back extension to release pressure in each spinal vertebra. ./.

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